

Gauze and Cotton tissue

一、 Product name, model and specification

Product name: Gauze and Cotton tissue

Model: dry type; Wet mold

Specifications: 5×6cm, 6×6cm, 12×15cm, 12×18cm, 13×15cm, 15×18cm, 15×20cm, 16×20cm, 18×18cm, 18×20cm, 20×20cm, 20×25cm, 25×25cm

二、 the record person/production enterprise name, address/production address, contact information and after-sales service unit

1、 Name: Hunan Zhenheyikang Medical Instrument Co.,LTD

2、 Address: No.6 Building Jingxiang Energy, No.55 Xiaguang east road, Gaoxin district Xiangtan, Hunan

3、 Production address:NO.6 Building, NO.8 Building,5th floor Jingxiang Energy, No.55 Xiaguang east road,Gaoxin district Xiangtan, Hunan

4、 Contact number: +86 17711627084

5、 After-sales service unit: Hunan Zhenheyikang Medical Instrument Co.,LTD

三、 Production Registration Certificate No: Xiangtan Medical Products Administration Production Registration NO 20140004

四、 Product Registration Certificate No、 Product Technical Requirements No: XTFDA Certified 20230024

五. Product description, intended use

Use water-absorbing materials. Contains no disinfectant. Non-sterile for single use.

It is used to clean skin and wound surface.

六, use method: open the bag can be used.

七. Contraindications, precautions, warnings and tips

Product contraindications: Do not use if you are allergic to the product material

Avoid humidity and direct sunlight;

Damaged packaging prohibited use, expired products prohibited use;

Disposable use, please destroy after use;

八. Interpretation of graphics, symbols, abbreviations and other contents used in medical device labels



一次性使用



参考使用说明书



包装破损切勿使用

九、 Installation and use instructions or illustrations

Select the appropriate model and specification to use directly。

十、 Product maintenance and maintenance methods, special storage conditions, methods

The product should be stored in a non-corrosive gas and well-ventilated room, ordinary transport, in the process of transportation, storage should be kept clean and pollution-free, strictly prohibited rain and snow immersion。

十一、 Date of manufacture, expiration date or expiration date

Valid for three years, production batch number, production date, expiration date refer to the small package or qualification label.

十二、 Accessories list, including accessories, accessories, wear and tear replacement cycle and replacement method description

This product has no accessories

十三、 The date and version of the preparation or revision of the specification
This specification was revised in March 2023, version 2023.